Dantec EMG Systems Somatics ECT Systems Mag Venture TMS Systems Neuro Accessories



SONORAY Pty Ltd Unit 14, 12a Loyalty Rd NORTH ROCKS

Phone: +61 2 9630 9638 **Fax:** +61 2 9630 5491

NSW 2151 AUSTRALIA

A.B.N. 91 057 148 973

All Mail to: PO Box 4823 NORTH ROCKS NSW 2151 AUSTRALIA

Email: sonoray@ozemail.com.au www.sonoray.com




Dantec DCN New Improved Design

Disposable Concentric Needle Electrodes for EMG

Premium performance with enhanced comfort

- Silicone coating and ultra-sharp, 3-bevel tip design increase patient comfort by reducing friction and insertion force
- Stainless steel cannula provides strength for efficient penetration and flexibility for precise adjustments
- Each needle has an ergonomically designed, color-coded hub with a recording area direction indicator
- Precisely centered, insulated tungsten core provides accurate signal detection
- Rigorous quality control assures consistent sharpness and high concentricity for clinical reliability
- Low impedance levels for accurate recordings
- Pre-sterilized





Tungsten core assures reliable and accurate signal detection

Recording area 0.02 or 0.07 mm² depending on needle size